H50-007-2020-11

CE / CE₀₄₈₂

Not for US-Distribution

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Warnings:

- The reuse of LINK single-use products is not permissible.
- The package inserts provided with the devices do not contain all of the information necessary for the selection and application of the devices. For proper handling, refer to other product-related instructions, such as those for the surgical technique associated with the relevant system and, if necessary, any special handling recommendations.
- Products made of plastic (e.g. polyamide (PA), polyethylene (PE), polyoxymethylene (POM), ultra high molecular weight polyethylene (UHMWPE)) may not be localizable using external imaging procedures.
- For the processing of LINK instruments, it is presumed that the personnel have technical knowledge level I (Germany) and in other countries technical knowledge and expertise.
- Medical devices that are sent in for servicing must be processed beforehand in such a way that they cannot constitute a hazard to third parties.

1. First use

Non-sterile instruments must be cleaned, disinfected and sterilized according to these instructions before first use. Before reprocessing an instrument, the persons in charge of the reprocessing unit for medical devices (AEMP) must compare the compatibility of their established reprocessing procedures with those described in this reprocessing instruction. All users in the OR must familiarize themselves with the function of an instrument prior to use.

2. Preparation at the site of use and transport to the reprocessing area

The instruments must be cleaned, disinfected, inspected and sterilized before their first use and after each use. For this purpose, they are to be disassembled into their individual parts as far as is intended. Coarse contamination of the instruments must be removed with a lint-free cloth immediately after use to prevent body fluids, tissue residues and other residues from drying on the instruments.

Suitable transport containers must be used for the safe transport of the instruments to the reprocessing area to protect the medical device, the environment and the medical staff. Careful and proper handling of the instruments must be ensured during transport. Particular care must be taken with sensitive instruments and instruments with cutting edges to prevent damage to the product.

The instruments must be disposed off dry and care must be taken not to exceed a period of six hours between use and reprocessing.

3. Inspection and manual pre-cleaning

In individual cases, refer to the additional, instrument-specific LINK cleaning and maintenance instructions.

All instruments must be pre-cleaned manually as a matter of principle; deviations from this are only permitted if no residues (e.g. blood residues or bone chips) are detected visually. For this purpose the instruments are removed from the trays, disassembled and rinsed under running water (temperature < 30°C, drinking water quality). A soft and suitable brush is used to clean surfaces, crevices, lumens and other difficult-to-reach areas until the visible residues have been removed.

All instruments must be cleaned ultrasonically as a matter of principle; this may only be deviated from if one can rule out that residues may still be present in areas that cannot be checked visually. The ultrasonic bath must be filled according to the manufacturer's instructions and a cleaning agent (mild alkaline detergent; suitable for ultrasonic cleaning) added to the water. Follow the manufacturer's instructions when determining the cleaning solution concentration and the water temperature.

The instruments are placed on trays and must be completely immersed in the cleaning solution.

Recommended cleaning conditions:

- Ultrasonic frequency: 25-40 kHz

After ultrasonic cleaning, the instruments must be rinsed with potable water for at least one minute to remove any detergent residues. The instruments must be inspected visually for cleanliness. The above described procedure must be repeated until no more visible contamination is present.

Deviations from the procedure prior to machine cleaning and thermal disinfection (e.g. omission or addition of individual manual reprocessing steps) for individual products are possible if the effectiveness of cleaning has been proven at the operator's site during validation of the reprocessing processes.

4. Machine cleaning and disinfection

We recommend the mechanical method for safe cleaning and disinfection. The processes of the WDs must be validated at the operator's.

For machine cleaning and disinfection, it is necessary to remove the instruments from the trays and place them in suitable trays disassembled and/or opened. The instruments must not protrude beyond the tray. Rotation of the spray arm must not be obstructed or blocked. Lightweight and/or small components must be secured on the trays such that they cannot fall out. Care must be taken to ensure that the instruments are rinsed sufficiently and that no areas are missed. Place medical devices on the load carrier of the WD such that all internal and external surfaces of the medical devices are cleaned and disinfected. Hollow bodied instruments must be positioned at a suitable point on the load carrier so that complete cleaning and disinfection of the internal surfaces of the medical devices is ensured. The WD operating manual contains information on the manner of loading and the correct loading procedure (ISO 15883-1, Section 8.3, b). Please also observe the information provided by the WD manufacturer.

Equipment: washer-disinfector according to ISO 15883 Part 1 and Part 2 with thermal disinfection

Machine cleaning and disinfection must be performed in a qualified WD (according to ISO 15883-1 and -2). We recommend the following program parameters for automatic cleaning and disinfection.

Program step	Water	Temperature	Time	Chemical
Pre-cleaning	Cold (at least potable water quality)	---	5 min	---
Cleaning	DIW	55 °C	10 min	Mildly alkaline
if necessary neutralization according to the manufacturer's instructions of the cleaning chemicals				
Post-rinsing	DIW	According to the manufacturer's specifications of the cleaning chemicals	2 x 1 min	---
Disinfection	---	A ₀ value 3000		
		90 °C	5 min	---
Drying	---	130°C	15 min	---

5. Inspection, maintenance and care

After each reprocessing procedure, the instruments must be checked for their operational and instrument-specific functionality. In particular, they must be free of visible residues and/or contamination (the use of an illuminated magnifying glass is recommended). Especially lumens and difficult to access areas must be checked. Principally, functions such as measuring functions, compatibilities, cutting edges, tips, connections, locks, notches and moving components are to be checked. Plastic components and instruments must be checked for wear and tear due to ageing as well as for cracks, embrittlement or chipping. Rotary instruments must be checked for deformation. The disassembled instruments must be assembled and checked for functional testing.

Damaged instruments must be replaced with new or functional instruments.

Perform an inspection or functional test of the instruments after each reprocessing procedure. To reduce friction and wear after cleaning and thermal disinfection, hinge joints, screw threads and other moving parts must be lubricated with an oil suitable for the sterilization process used (spray can, oil pen or dropper bottle) (one to two applications). In individual cases, also refer to the additional, instrument-specific LINK cleaning and maintenance instructions.

Note: the care product must be based on paraffin/white oil, must be biocompatible, steam sterilizable and steam-permeable according to the valid European or United States Pharmacopoeias.

The service life of instruments depends on the material, design, application and reprocessing. Therefore, do not use damaged instruments, instruments with visible surface changes or instruments whose labeling is no longer clearly recognizable.

Note: do not carry out any repairs yourself. Service and repairs may only be carried out by suitably qualified personnel of the manufacturer.

6. Packaging

Cleaned and disinfected instruments, in assembled condition, must be packaged with a validated sterile barrier system for sterilization. Ensure that the packaging is compatible with the sterilization process.

7. Sterilization

LINK instruments are designed for sterilization with moist heat (steam sterilizer according to EN 285) using the fractionated pre-vacuum method. The sterilization process according to ISO 17665 must be validated by the user.

We recommend the following program parameters for the sterilization process:

	Europe
Procedure	Fractionated pre-vacuum process
Temperature	134 °C
Hold time	At least 3 min
Drying time	Recommended: 20 min Method: drying vacuum

Further information on the proper reprocessing of medical devices can also be obtained from the „Arbeitskreis Instrumentenaufbereitung“ [Working Committee Instrument Reprocessing] at www.a-k-i.org.

8. Storage

Surgical instruments must always be handled with care. This particularly applies during transport, cleaning, care, sterilization, and storage. Following sterilization, the instruments must be stored in a dust-free and dry place. The storage period must be determined in detail with the person responsible for hygiene at the operator's. Avoid direct sunlight. Improper handling or care as well as misuse can lead to premature wear or damage.

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